

Neutral pH Detergent

For All Brands of Automatic Washers - Low Foam

Product Features

Cleans a wide range of metal surfaces - Safe and effective for stainless steel, aluminum, tungsten carbide, copper and brass instruments as well as mixed metals and all other materials.

Free of metal-attacking agents - does not contain scale-forming phosphates, chlorides or hydroxides.

Counteracts hard water - chelating additives help protect steel by preventing scale build up.

Penetrates bio-burden - frees instruments of stubborn organic encrustation.

Versatile - formulated for all brands of automatic washers/disinfectors.

Directions for Use

CST-400 can be used manually and automatically at 5-20 ml per gallon (1-5ml per liter). For most cleaning applications, 5 ml per gallon is sufficient to ensure effective cleaning. Occasionally, more **CST-400** may be required depending upon water quality and bio-burden load. Make sure that machine goes through both wash and rinse cycles. Post inspect all instruments. If soil or residue is present, reprocess using a higher level of **CST-400**. For manual use, discard spent wash solution on a daily basis.

Specifications

Chemical Description: Concentrated liquid blend of non-ionic surfactants and chelating agents

Physical Appearance: Clear, low viscosity liquid, easily drawn through automatic metering systems

pH: 6-8 at recommended dilution ratio 6.5-7.8 as is

Specific Gravity: 1.03 g/ml

Order Information

CST-400-1 - 1 gallon plastic bottle packed four per case

CST-400-5 - 5 gallon plastic jug with handle

CST-400-15 - 15 gallon plastic drum

CST-400-30 - 30 gallon plastic drum

Cautions/Warnings

Read Material Safety Data Sheet (MSDS) before using this product. KEEP OUT OF REACH OF CHILDREN. As a safety precaution, protect skin and eyes by wearing NIOSH approved safety goggles and gloves.



COMPLETE SOLUTIONS TECHNOLOGIES
CST

Complete Solutions Technologies, LLC
302 N. Washington St., Suite 202 E
Moorestown, NJ 08057

www.cstllc.com
Toll-Free 888-987-1987
Tel 856-234-7392
Fax 856-234-7918