

Neutralizer

NEUTRALIZES ALKALINE DETERGENT IN AUTOMATIC WASHERS

CST-200 - when used in the rinse cycle of automatic washer/decontaminators, ensures complete neutralization of residual alkaline detergent carry over from the cleaning cycle.

Product Features

Provides additional instrument protection - acidic rinse promotes natural passivation process.

Fights corrosion - anti-corrosive agents inhibit corrosion which help extend instrument life.

Hard water management - chelating additives helps protect instruments against spotting & staining.

Multi-functional - can be use as acidic cleaner.

Directions for Use

Program machine to deliver one tenth of a fluid ounce per gallon water (rinse cycle).

Technical Data

Chemical Description: Natural, organic acid liquid concentrate with chelating additives and corrosion inhibitors

Physical Appearance: Clear water white liquid with a faint odor

pH: 1.70

Specific Gravity: 1.2

Order Information

- CST-200-5** - 5 gallon plastic jug
- CST-200-15** - 15 gallon plastic drum
- CST-200-30** - 30 gallon plastic drum
- CST-200-55** - 55 gallon plastic drum

Caution/Warnings

Care should be taken to minimize exposure. As a safety precaution, wear NIOSH approved safety gloves and goggles. Refer to MSDS for additional safety information.

