

Lubricant

FOR USE IN AUTOMATIC WASHERS

CST-300 - when used in the lubrication cycle of automatic washer/decontaminators, helps protect surgical instruments from corrosion and staining associated with the harsh elements of steam sterilization.

Product Features

Gentle formula - formulated with gentle, non-aggressive ingredients.

Extends instrument life - leaves a protective microfilm which prolongs instrument life.

Frees up box locks - lubricates moving parts of surgical instruments.

Economical - high level of concentration is cost effective.

Directions for Use

Program machine to deliver Only 1/10 of one ounce (3 ml) per gallon water.

Technical Data

Chemical Description: Water emulsible, non-ionic blend of corrosion inhibitors and gentle lubricants

Physical Appearance: Clear, light yellow liquid with a slight scent

pH: 6.8 - 7.5 at recommended dilution ratio

Specific Gravity: 0.89

Order Information

CST-300-5 - 5 gallon plastic jug
CST-300-15 - 15 gallon plastic drum

Caution/Warnings

Care should be taken to minimize exposure. As a safety precaution, wear NIOSH approved safety gloves and goggles. Refer to MSDS for additional safety information.

KEEP OUT OF REACH OF CHILDREN!



COMPLETE SOLUTIONS TECHNOLOGIES
CST

Complete Solutions Technologies, LLC
 302 N. Washington St., Suite 202 E
 Moorestown, NJ 08057

www.cstllc.com
1-888-987-1987
Tel 856-234-7392
Fax 856-234-7918