

INSTRUMENT SPRAY LUBRICANT

For Convenient Protection of Surgical Instruments

CST-302S is a convenient, ready-to-use spray lubricant, specially formulated to protect and extend the working life of surgical instruments.

Product Features

Helps prevent corrosion – anti-corrosive agents fight corrosion.

Gentle formula – formulated with gentle, non-aggressive ingredients.

Extends instrument life – leaves a protective micro-film which prolongs instrument life.

Frees up box locks – lubricates moving parts of surgical instruments.

Convenient to use – sprayer feature eliminates time-consuming milk

Directions for Use

Spray **CST-302S Instrument Spray Lubricant** onto clean stainless steel surgical instruments. Allow to drip dry, then proceed to autoclave.

Technical Data

Chemical Description:	Water emulsible, non-ionic lubricant emulsion milk
Physical Appearance:	Liquid milk appearance with pleasant scent
pH:	6.8-7.2
Specific Gravity	0.99

Order Information

CST-302S - 8	8 fl.oz. HDPE sprayer bottles packed 12 per case
CST-302S - 32	32 fl.oz. HDPE sprayer bottles packed 12 per case

Cautions/Warnings

Instruments must be sterilized following treatment with **CST-302S**; follow your standard sterilization procedures.



COMPLETE SOLUTIONS TECHNOLOGIES
CST

Complete Solutions Technologies, LLC
302 N. Washington St., Suite 202 E
Moorestown, NJ 08057

www.cstllc.com
Toll-Free 888-987-1987
Tel 856-234-7392
Fax 856-234-7918