

Two-in-One Instrument Cleaner and Lubricant

Instrument Detergent/Lubricant for Multiple Washer Types

Product Features

Cleans a wide range of metal surfaces - safe and effective for stainless steel, aluminum, tungsten carbide, copper and brass instruments as well as mixed metals and many other materials.

Free from metal-attacking agents - does not contain scale-forming phosphates, chlorides, hydroxides; formula is pH balanced and features chelating additives to counteract hard water.

Highly effective bio-burden penetrant - frees instruments of stubborn organic encrustation.

Versatile - formulated for manual, ultrasonic including laboratory glassware cleaning.

Lubricant - lubricates and keeps box locks open.

Free rinsing - formula has the ability to permit rinsing of instruments to a state of total freedom from any residue.

Directions for Use

CST-304 can be used manually or automatically at 5-20 ml. per gallon (1/6 to 2/3 fl.oz. per gallon. This quantity is sufficient to ensure effective cleaning for most applications, though on occasion more **CST-304** may be required. Ensure machine completes both cycles, post inspect, and if residue remains, reprocess using a higher level of **CST-304**.

Specifications

Chemical Description: Liquid surfactant blend with chelating agents and mild scent

Physical Appearance: Viscous yellow liquid, easily dispensed with manual or automatic pump

pH: 7.0 - 8.0 at recommended dilution ratio

Specific Gravity: 1.03 g/ml

Order Information

- CST-304-1** - 4 x 1 gallon bottles per case
- CST-304-5** - 5 gallon plastic jug with handle
- CST-304-15** - 15 gallon plastic drum
- CST-304-30** - 30 gallon plastic drum
- CST-304-55** - 55 gallon plastic drum

Cautions/Warnings

Read Material Safety Data Sheet (MSDS) before using this product. **KEEP OUT OF REACH OF CHILDREN!**

As a safety precaution, protect skin and eyes by wearing NIOSH approved safety goggles and gloves.



COMPLETE SOLUTIONS TECHNOLOGIES
CST

Complete Solutions Technologies, LLC
302 N. Washington St., Suite 202 E
Moorestown, NJ 08057

www.cstllc.com
Toll-Free 888-987-1987
Tel 856-234-7392
Fax 856-234-7918