

Extreme Alkaline Detergent

For Automatic Washers/Decontaminators

Product Features

Optimal buffered alkaline system - formulated with powerful buffers that maintain a high pH level, even in severe hard water conditions.

Agressively removes blood and bio-burden - alkalinity is excellent for maintaining stainless steel instruments.

Compatible with all automatic washers/decontaminators - easily integrated with all autowasher injection systems configured for alkaline cleaning.

Economical - high-end cleaning results at competitive cost.

Prevents hard water deposits - superior preventative action is achieved by sequestering and complexing hard water ions.

Directions for Use

Program the machine to deliver 4 ml/liter (0,0135 fl.oz./gal.) of **CST-101**. On occasion, the quantity of **CST-101** may be adjusted contingent upon bio-burden load and water quality. A pH of 11.2 - 11.8 is a reliable indicator of the correct dosage. Use **CST-200** at 2ml/liter to neutralize residual alkaline detergent.

Specifications

pH:	11.6 - 11.9
Specific Gravity:	1.5 g/ml
Conductivity:	2.09 ms/cm

Order Information

CST-101-5	-	5 gallon plastic jug with handle
CST-101-15	-	15 gallon plastic drum
CST-101-30	-	30 gallon plastic drum
CST-101-55	-	55 gallon plastic drum

Cautions/Warnings

CST-101 is highly alkaline and can cause burning of the skin and eyes. As a safety precaution, protect skin and eyes by wearing NIOSH approved safety goggles and gloves. Refer to MSDS for additional safety information.
KEEP OUT OF REACH OF CHILDREN!



COMPLETE SOLUTIONS TECHNOLOGIES
CST

Complete Solutions Technologies, LLC
302 N. Washington St., Suite 202 E
Moorestown, NJ 08057

www.cstllc.com
Toll-Free 888-987-1987
Tel 856-234-7392
Fax 856-234-7918